

CHRONIC WOUND

CONSUMABLES SCHEME (CWCS)

MEDSTOCK CWCS-LISTED PRODUCTS

All dressings below are fully subsidised through the CWCS. Designed for a range of wound types, our dressings offer outstanding value and clinical performance.

No-Cost Wound Dressings for Eligible Patients Available through the Chronic Wound Consumables Scheme (CWCS)

Visit: www.health.gov.au/our-work/cwcs

Not registered to prescribe? Email: cwcs-training@monash.edu

PRIMARY DRESSINGS

FOAM NON-ADHESIVE

Suitable for pressure injuries, vascular ulcers, partial thickness burns and mixed aetiology ulcer

- + Soft, absorbent, and can be cut to size
- + Contains CMC granules that expand on contact with moisture to fill dead space and manage moderate to high exudate levels
- + Vertically wicks to protect against maceration



Product Code: MS0505F
Product Size: 5cm x 5cm
Pcs/Box: 10



Product Code: MS1010F
Product Size: 10cm x 10cm
Pcs/Box: 10



SILICONE FOAM LITE WITH BORDER

Suitable for skin tears, low exuding wounds, ulcers, pressure injuries, blisters, abrasions and fragile skin

- + Thin profile 4-layered gentle adhesive foam dressing
- + Promotes a moist wound healing environment




Product Code: MS0505LIT
Product Size: 5cm x 5cm
Pcs/Box: 10



Product Code: MS7575LIT
Product Size: 7.5cm x 7.5cm
Pcs/Box: 10



Product Code: MS1010LIT
Product Size: 10cm x 10cm
Pcs/Box: 10



SILICONE FOAM WITH BORDER

Suitable for skin tears, vascular ulcers, diabetic ulcers, pressure injuries, and fragile skin

- + Protective 5-layered gentle adhesive foam dressing
- + Ideal for moderate to highly exuding wounds
- + Promotes a moist wound healing environment



Product Code: MS7575SF
Product Size: 7.5cm x 7.5cm
Pcs/Box: 10



Product Code: MS1010SF
Product Size: 10cm x 10cm
Pcs/Box: 10



ALGINATE

Suitable for vascular ulcers, sinuses, open cavities and skin donor sites

- ✦ Made from high-mannuronic acid (High M) 100% calcium alginate fibres
- ✦ Absorbs moderate to high levels of exudate
- ✦ Forms gel to promote a moist environment
- ✦ May promote haemostasis in minor bleeding



Product Code: MS51010A
Product Size: 10cm x 10cm
Pcs/Box: 5



Product Code: MS1010A
Product Size: 10cm x 10cm
Pcs/Box: 10



Product Code: MS5230SA
Product Size: 2cm x 30cm
Pcs/Box: 5



Product Code: MS0230A
Product Size: 2cm x 30cm
Pcs/Box: 10

Up to
7
DAYS WEAR TIME

HIGH
EXUDATE


NON-ADHESIVE

PRIMARY OR SECONDARY DRESSINGS

SUPER ABSORBER

Suitable for pressure injuries, vascular ulcers, skin grafts and burns

- ✦ For use with high to very highly exuding wounds
- ✦ Requires fixation to secure
- ✦ Effective under class I and class II compression
- ✦ Ideal for use over cavities or sinuses
- ✦ Pair with alginate rope for enhanced cavity management
- ✦ Absorbs up to 40 times its dry weight



Product Code: MS1010SA
Product Size: 10cm x 10cm
Pcs/Box: 10



Product Code: MS1015SA
Product Size: 10cm x 15cm
Pcs/Box: 10

Up to
7
DAYS WEAR TIME

V.HIGH
EXUDATE

NON-ADHESIVE

MEDSTOCK CWCS-LISTED PRODUCTS

The Chronic Wound Consumables Scheme (CWCS) fully subsidises the cost of wound consumable products for people who:

- ✦ Have diabetes AND a chronic wound (present for >6 weeks), AND
- ✦ Are 65 years and over, OR 50 years and over for First Nations people, AND
- ✦ Hold a valid Medicare card.

This subsidy cannot be used in conjunction with:

- ✦ The National Disability Insurance Scheme (NDIS)
- ✦ The Department of Veterans' Affairs rehabilitation and treatment programs, including Repatriation Pharmaceutical Benefits Scheme (RPBS) dressings
- ✦ Australian Government-funded residential aged care homes, Home Care Packages, or the Community Home Support Program
- ✦ Public or community wound clinics run by or on behalf of state and territory governments that fund wound care products

Healthcare professionals must be trained and registered with CWCS to prescribe consumables.

WHY CHOOSE MEDSTOCK?



Proudly Australian-owned and operated



ISO 9001 and 13485



Unique colour-coded packaging to reduce error



Extensive clinical support and education



Competitive pricing across all product lines